



TY001-Cuspidor

01TY01-00A

No.	Material No.	Name	Quantity	Remark
1	8016747A	φ5 Clamp	2	
2	8019006A	TY-1 Connecting Pin of Rotating Shaft	1	
3	8018997B	TY-1 Rotating Shaft	1	
4	8026038A	TY-1 Connecting Plate	1	
5	8025894A	M6 * 10 stainless steel slotted flat head screws	2	
6	8020173A	M5*16 Stainless Steel Socket Head SEM Screw	3	
7	8019018B	TY-1 Connecting Hose of Funnel	1	
8	8026044A	φ6.9 Copper Connector	1	
9	8019591A	ST2.2*10 Stainless Steel Self-tapping Philips Screw	6	
10	8019393A	Rubber Plug	6	
11	8022600A	TY-1 Cover of Assistant Arm	1	
12	8010858A	3/8 inch yellow tube	1	
13	8026062A	TY-1 under Shell	1	
14	8018229A	M6*30 Stainless Steel Socket Head SEM Screw	4	
15	8024361A	Air Pilot Valve 2	2	
16	8011032B	φ4.7 Copper Connector	4	
17	8026050A	Hex Three-Way Valve	2	
18	8012199B	φ2.4 Copper Connector	4	
19	8025897A	M10*1 Stainless Steel Nut	2	
20	8019450B	MNV-3 Pin Valve	2	
21	8026026A	Valve body fixing plate	1	
22	8003751A	M4*10 Stainless Steel Countersunk Head Philips Screw	2	
23	8019453A	2 Way Switch	2	
24	8026065A	TY-1 Upper Shell	1	
25	8009280A	M14*2 Nut	2	
26	8022492A	M14 Double Self-locking Washer	1	
27	8014227A	Seal Washer 15*23*3	1	
28	8011032B	φ4.7 Copper Connector	1	
29	8013750B	Water Hose Holder	1	
30	8014224A	O Ring of φ 11.8* φ 1.8	4	
31	8013747A	Water Supply Hose	1	
32	8003766A	ST4*16 Stainless Steel Pan Head Self-tapping Philips Screw	2	
33	8011032B	φ4.7Copper Connector	1	
34	8008707A	M4*6 Stainless Steel Socket Head Set Screw	1	
35	8026083A	Fixed Plate for Flush Hose	1	
36	8026068A	O Ring of φ5*φ1.5	2	
37	8026029B	TY-1Flush Hose	1	
38	8003817A	CCeramic Cuspidor Cover	1	
39	8008710A	Cuspidor Filter (45mm)	1	
40	8000916A	AJ15 Cuspidor Funnel	1	
41	8009061A	Cuspidor Rubber Mat	1	
42	8019048A	Circular Ceramic Cuspidor	1	
43	8007372C	Nut of AJ15 Cuspidor Funnel	1	
44	8007858A	φ54.5*φ3.55 O Ring	2	
45	8019015A	TY-1 Funnel	1	
46	8018994E	TY-1 Rotating Arm	1	
47	8001702A	M5*16 Stainless Steel Pan Head Philips Screw	4	
48	8019012A	TY-1 Cover of Rotating Arm	1	