



Ceiling light Amber LED Light Parts list

02KD02-00A

No.	Part number	Name	Number	Remark
1	8022876A	M5*8 Stainless steel screw	4	
2	8013801A	M6*8 set bolt (black)	1	
3	8001684A	M4*16 Stainless steel screw	1	
4	8021175A	M5*6 Stainless steel screw	1	
5	8017785B	ceiling light	1	
6	8022522A	copper sleeve	1	
7	8017785B	HK3520needle bearing	1	
8	8022522A	M6*8 set bolt (black)	1	
9	8017779C	light arm connecting shaft	1	
10	8022528A	DL-2 Rotation limit	1	
11	8022876A	M5x8 Stainless steel screw	4	
12	8017791B	adjustable light arm II-balance arm (1.72Kg)	1	
13	8023527A	DL-2 ceiling light rotation arm	1	
14	8011491A	glass cuspidor silicone decoration cover	2	
15	8001738A	M6*16 Stainless steel screw	2	
16	8001564A	φ6 spring pad	2	
17	8023530A	ceiling light balance arm fixed seat	1	
18	8023533A	ceiling light balance arm fixed seat seal board	1	
19	8020617A	Amber LED dental light (ADS)	1	
20	8014164A	M6*5 alloy steel set screw	1	