


# Classic Instrument Tray Breakdown

02GJP02-01A



## Classic Instrument Tray Breakdown

02GJP02-01A

No.	Material No.	Name	Quantity	Remark
1	8014998A	M5*12 Stainless Steel Truss Head Philips SEM Screw	1	
2	8024655A	SMDE-1 Bottle Fixed Plate	1	
3	8024754B	2-way Valve ( $\phi$ 2.4, $\phi$ 4.7)	1	
4	8018112A	Mini Pressure Regulating Valve	1	
5	8023635B	Circular Manometer (Metric System)	1	
6	8018115A	G1/8 Hex End Cap	1	
7	8008719E	Bottle Head (with installation hole)	1	
8	8013687A	Rubber Mat ( $\phi$ 28* $\phi$ 18*t2)	2	
9	8018931B	Fixed Axis of Rotatory Bottle	2	
10	8018904B	Connecting Shafts of Rotatory Bottle	1	
11	8020602A	O Ring of $\phi$ 20* $\phi$ 2.4	1	
12	8018907B	Head of Rotatory Bottle	1	
13	8026344A	Ajax Water Bottle 1.5L	1	