


Whale LED Light with arm and bushing Parts list

100KD02-00A



No.	Material No.	Name	Quantity	Remark
1	8013552C	Whale dental light	1	
2	8013558E	light arm	1	
3	8020503A	Power Cord	1	
4	8020560B	Lower Shell of Power Box	1	
5	8015025A	Fuseholder R3-53A	1	
6	8015019A	6GFU-F6A250V Fuse	1	
7	8015040A	2 Ways Block (DMTF250)	1	
8	8020563A	30VA Transformer	1	
9	8020557A	Upper Shell of Power Box	1	
10	8016762A	M4*6 Stainless Steel Truss Head Philips Screw	2	