



## 201 rear cuspidor bom diagram15

201ZS01-00A

No.	Part number	Number	Material	Remark
1	8026017A	rear cuspidor assistant arm	1	
2	8011584A	1/4 inch pipe quick connector	2	
3	8018127A	regulating valve	1	
4	8012196A	small nipple water regulating valve (φ2.4-φ2.4)	2	
5	8026035A	external water-air bracket shell	1	
6	8019912A	6#-32*3/8 Stainless steel screw	2	
7	8016660A	M6*20 Stainless steel screw	2	
8	8026032A	external water-air bracket shell	1	
9	8022177A	M4*12 Stainless steel screw	2	
10	8026020A	rear cuspidor arm shaft sleeve	1	
11	8017272B	filter (Large straight out)	1	
12	8018229A	M6*30 Stainless steel screw	8	
13	8018943E	rear tank bracket	1	
14	8012721A	AJ16 front rotate arm limit sleeve	4	
15	8001738A	M6*16 Stainless steel screw	4	
16	8022426A	304 Stainless steel with glue hoop	1	
17	8009403A	M5*10 Stainless steel screw	2	
18	8021163A	M10*30 Stainless steel screw	4	
19	8026053A	rear cuspidor connection board	1	
20	8011842A	Φ50 down connection shaft	2	
21	8011830A	AXK5070 surface bearing	4	
22	8011842A	φ50 connection shaft cover	2	
23	8021163A	M10*30 Stainless steel screw	2	
24	8025891A	Φ52X45 grey corrugated pipe (length 80CM)	1	