


16 US Junction Box Breakdown

300DX01-00A



No.	Material No.	Name	Quantity	Remark
1	8017680A	AJ16 Junction Box Cover II	1	
2	8015112A	M4*8 Stainless Steel Truss Head Philips Screw	4	
3	8018109A	Air Filter Pressure Reducing Valve	1	
4	8018106A	Water Filter Pressure Reducing Valve	1	
5	8001252A	Six-way Connector	1	
6	8021664A	Three-way Connector (φ2.4)	1	
7	8010840D	AJ16 Junction Box Base Plate	1	
8	8012160A	AJ16 Fixed Holder of Front Shell	2	