



AJ12



# AJ12 Chair

**Part number: A091502**

Electromechanical chair with 3 programmable positions.

5 year warranty.



# AJ12 Chair

- AJ12 Electromechanical chair with 3 programmable positions,
- Articulating headrest,
- Swing out armrest,
- Steel backrest stand,
- Steel base plate,
- Chair safety switch,
- Wide backrest, Upgraded fiber leather,
- One touchpad attached on armrest, Foot control

## Specifications:





Environmental requirements — Ambient temperature  $5^{\circ}\text{C}\sim 40^{\circ}\text{C}$ ;  
 Humidity  $\leq 80\%$ ;  
 Atmospheric pressure 8.7 psi ~ 15.4psi;

Electrical requirements — AC  $110\text{V} \pm 22\text{V}$ ;  
 50HZ, 60HZ;

Motor run time — Intermittent operation 30s on / 5min Off;

Weight — 320lbs;

Max patient load — 400lbs;



# AJ12 Classic 100

Part number: A9121001



## Standard Features:

AJ12 Electromechanical chair with 3 programmable positions,  
Articulating headrest, Swing out armrest,  
Steel backrest stand, Steel base plate, Chair safety switch,  
Wide backrest, Upgraded fiber leather,  
One touchpad, Foot control  
Classic dentist's control head with 3 auto H.P. control,  
Individual water coolant controls,  
Master handpiece air coolant adjustment,  
Asepsis handpiece tubing,  
Master on/off switch, 1 doctor syringe, Wet/Dry foot control,  
Stainless steel tray with non-slip pad, Flex arm air brake switch

Classic 100 Post mounted utility (PMU) with steel side box,  
City water/bottled system, Water outlet and air QD,  
Solids collector,  
2 HVE, 1SE, 1 syringe on telescoping assistant's control,  
Junction box with air/water filters, regulators, master shut-off valves,  
5 ft. umbilical  
Post mount LED light with sensor control, 3 level intensity adjustment

# AJ12 Classic 101

Part number: A9121011



## Standard Features:

AJ12 Electromechanical chair with 3 programmable positions,  
Articulating headrest, Swing out armrest,  
Steel backrest stand, Steel base plate, Chair safety switch,  
Wide backrest, Upgraded fiber leather,  
One touchpad, Foot control  
Classic dentist's control head with 3 auto H.P. control,  
Individual water coolant controls,  
Master handpiece air coolant adjustment,  
Asepsis handpiece tubing,  
Master on/off switch, 1 doctor syringe, Remote toggles for cup fill  
and bowl rinse, Wet/Dry foot control,  
Stainless steel tray with non-slip pad, Flex arm air brake switch,  
Porcelain cuspidor bowl with strainer, Cup fill & Bowl rinse toggles.

Classic 101 Post mounted utility (PMU) with steel side box,  
City water/bottled system, Water outlet and air QD,  
Solids collector,  
2 HVE, 1SE, 1 syringe on telescoping assistant's control,  
Junction box with air/water filters, regulators, master shut-off valves,  
5 ft. umbilical  
Post mount LED light with sensor control, 3 level intensity adjustment

# AJ12 Classic 200

Part number: A9122001



## Standard Features:

AJ12 EL chair with 3 programmable positions,  
Articulating headrest, Swing up/down aluminium armrest,  
Narrow backrest, upgraded fiber leather, Cast iron chair frame,  
Steel backrest stand, Steel base plate,  
60 degrees seat rotation, chair safety switch,  
One touchpad, Foot control  
Swing Classic dentist's control head with 3 auto H.P. control,  
Individual water coolant controls,  
Master handpiece air coolant adjustment,  
Asepsis handpiece tubing,  
Master on/off switch, 1 doctor syringe,  
Wet/Dry foot control, Stainless steel tray with non-slip pad,  
City/bottled water system, 5 ft umbilical,  
Junction box with air/water filters, regulators, master shut-off valves,

Chair mounted swing rear assistant vacumm package,  
Water outlet and air QD on assistant tray,  
Solids collector,  
2 HVE, 1SE, 1 syringe  
Post mount LED light with sensor control, 3 level intensity  
adjustment



# AJ12 Classic 201

Part number: A9122011



## Standard Features:

AJ12 EL chair with 3 programmable positions,  
Articulating headrest, Swing up/down aluminium armrest,  
Narrow backrest, upgraded fiber leather, Cast iron chair frame,  
Steel backrest stand, Steel base plate,  
60 degrees seat rotation, chair safety switch,  
One touchpad, Foot control  
Swing Classic dentist's control head with 3 auto H.P. control,  
Individual water coolant controls,  
Master handpiece air coolant adjustment,  
Asepsis handpiece tubing,  
Master on/off switch, 1 doctor syringe,  
Wet/Dry foot control, Stainless steel tray with non-slip pad,  
City/bottled water system, 5 ft umbilical,  
Junction box with air/water filters, regulators, master shut-off valves,

Swing mount Cupspidor and assistant vacuum package,  
City water/bottled system, Water outlet and air QD,  
Solids collector,  
2 HVE, 1SE, 1 syringe on telescoping assistant's control,  
Junction box with air/water filters, regulators, master shut-off valves,  
5 ft. umbilical  
Post mount LED light with sensor control, 3 level intensity adjustment


# Dental Stools



## D3 Doctor stool

Part number: A080103


Chrome plated pneumatic cylinder  
 Height adjustment,  
 Backrest forward backward adjustment,  
 Seat surface adjustment,  
 Contoured seat & back,  
 Fiber leather upholstery  
 Cast aluminum star,  
 Quiet twin wheel casters,  
 1 year warranty




## N4 Assistant stool

Part number: A080224

Round seat,  
 Fiber leather upholstery,  
 Chrome plated pneumatic cylinder  
 Adjustable chrome foot ring  
 Cast aluminum star  
 Ratcheting body support,  
 1 year warranty



Color sample  
 (Fiber leather)



Technical support: 800 488 9708

Add: 1590 S Milliken Ave., unit A, Ontario, CA 91761

Tel: (626) 620 0456

[www.adsdental.com](http://www.adsdental.com) [www.adsequip.com](http://www.adsequip.com)

[sales@adsequip.com](mailto:sales@adsequip.com)