



ADS Dental System



# AJ17 Pediatric Chair

# AJ17 Pediatric Chair

Part number: A091701



## Standard Features:

AJ17 Electromechanical Pediatric chair, telescoping backrest,  
 Little size for 3-10 year old kids,  
 3 programmable positions,  
 0 position and LP last position,  
 Articulating headrest,  
 Swing out armrest,


Steel backrest stand,  
 Steel base plate,  
 Chair safety switch,  
 Upgraded fiber leather,  
 One touchpad attached on armrest, Foot control

## Specifications:

Environmental requirements — Ambient temperature 5 °C~40 °C;  
 Humidity ≤80%;  
 Atmospheric pressure 8.7 psi~ 15.4psi;  
 Electrical requirements — AC110V ±22V;  
 50HZ, 60HZ;  
 Motor run time — Intermittent operation 30s On / 5min Off;

Weight — 265lbs;

Max patient load — 400lbs;



# AJ17 Classic 100

Part number: A9171001



## Standard Features:

AJ17 Electromechanical Pediatric chair with telescoping backrest,  
For 3-10 year old kids, 3 programmable positions,  
Articulating headrest, Swing out armrest, Steel backrest stand,  
Steel base plate, Chair safety switch, Upgraded fiber leather,  
One touchpad, Foot control  
Classic dentist's control head with 3 auto H.P. control,  
Individual water coolant controls,  
Master handpiece air coolant adjustment,  
Asepsis handpiece tubing,  
Master on/off switch, 1 doctor syringe,  
Remote toggles for cup fill and bowl rinse, Wet/Dry foot control,  
Stainless steel tray with non-slip pad, Flex arm air brake switch,  
Porcelain cuspidor bowl with strainer, Cup fill & Bowl rinse toggles.

Classic 100 Post mounted utility (PMU) with steel side box,  
City water/bottled system, Water outlet and air QD,  
Solids collector,  
2 HVE, 1SE, 1 syringe on telescoping assistant's control,  
Junction box with air/water filters, regulators, master shut- off valves,  
5 ft. umbilical,  
Post mount LED light with sensor control, 3 level intensity adjustment

# AJ17 Classic 101

Part number: A9171001



## Standard Features:

AJ17 Electromechanical Pediatric chair with telescoping backrest,  
For 3-10 year old kids, 3 programmable positions,  
Articulating headrest, Swing out armrest, Steel backrest stand,  
Steel base plate, Chair safety switch, Upgraded fiber leather,  
One touchpad, Foot control  
Classic dentist's control head with 3 auto H.P. control,  
Individual water coolant controls,  
Master handpiece air coolant adjustment,  
Asepsis handpiece tubing,  
Master on/off switch, 1 doctor syringe,  
Remote toggles for cup fill and bowl rinse, Wet/Dry foot control,  
Stainless steel tray with non-slip pad, Flex arm air brake switch,  
Porcelain cuspidor bowl with strainer, Cup fill & Bowl rinse toggles.

Classic 101 Post mounted utility (PMU) with steel side box,  
City water/ bottled system, Water outlet and air QD,  
Solids collector,  
2 HVE, 1SE, 1 syringe on telescoping assistant's control,  
Junction box with air/water filters, regulators, master shut- off valves,  
5 ft. umbilical  
Post mount LED light with sensor control, 3 level intensity adjustment

# AJ19 Pediatric Bench

Part number: A091901



## Standard Features:

AJ19 Electromechanical pediatric bench with up down function, Foot control articulating headrest  
 Articulating headrest, Sturdy cast iron structure & aluminum base,  
 Chair safety switch, Foot control .

## Specifications:





Environmental requirements — Ambient temperature  $5^{\circ}\text{C} \sim 40^{\circ}\text{C}$ ;  
 Humidity  $\leq 80\%$ ;  
 Atmospheric pressure  $8.7 \text{ psi} \sim 15.4 \text{ psi}$ ;

Electrical requirements — AC  $110\text{V} \pm 22\text{V}$ ;  
 50HZ, 60HZ;

Motor run time — Intermittent operation 30s On / 5min Off;

Weight — 331lbs;

Max patient load — 400lbs;



# AJ19 Classic 100

Part number: A9191001



## Standard Features:

AJ19 Electromechanical pediatric bench with up down function,

Articulating headrest,

Chair safety switch,

Foot control

Classic dentist's control head with 3 auto H.P. control,

Individual water coolant controls,

Master handpiece air coolant adjustment,

Asepsis handpiece tubing,

Master on/off switch, 1 doctor syringe,

Wet/Dry foot control,

Stainless steel tray with non-slip pad,

Flex arm air brake switch

Classic 100 Post mounted utility (PMU) with steel side box,

City water/bottled system, Water outlet and air QD,

Solids collector,

2 HVE, 1SE, 1 syringe on telescoping assistant's control,

Junction box with air/water filters, regulators, master shut- off valves,

5 ft. umbilical,

Post mount LED light with sensor control, 3 level intensity adjustment

# AJ19 Classic 101

Part number: A9191011



## Standard Features:

AJ19 Electromechanical pediatric bench with up down function,

Articulating headrest,

Chair safety switch,

Foot control

Classic dentist's control head with 3 auto H.P. control,

Individual water coolant controls,

Master handpiece air coolant adjustment,

Asepsis handpiece tubing,

Master on/off switch, 1 doctor syringe,

Wet/Dry foot control,

Stainless steel tray with non-slip pad,

Flex arm air brake switch

Classic 101 Post mounted utility (PMU) with steel side box,

City water/ bottled system, Water outlet and air QD,

Solids collector,

2 HVE, 1SE, 1 syringe on telescoping assistant's control,

Junction box with air/water filters, regulators, master shut- off valves,

5 ft. umbilical

Post mount LED light with sensor control, 3 level intensity adjustment


# Dental Stools



## D3 Doctor stool

Part number: A080103


Chrome plated pneumatic cylinder  
 Height adjustment,  
 Backrest forward backward adjustment,  
 Seat surface adjustment,  
 Contoured seat & back,  
 Fiber leather upholstery  
 Cast aluminum star,  
 Quiet twin wheel casters,  
 1 year warranty




## N4 Assistant stool

Part number: A080224

Round seat,  
 Fiber leather upholstery,  
 Chrome plated pneumatic cylinder  
 Adjustable chrome foot ring  
 Cast aluminum star  
 Ratcheting body support,  
 1 year warranty



Color sample  
 (Fiber leather)



Technical support: 800 488 9708

Add: 1590 S Milliken Ave., unit A, Ontario, CA 91761

Tel: (626) 620 0456

[www.adsdental.com](http://www.adsdental.com) [www.adsequip.com](http://www.adsequip.com)

[sales@adsequip.com](mailto:sales@adsequip.com)